
Scaut

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
User Manual

My intended use

I, Scaut, am a digital dermatoscope, developed to enable the visual examination and documentation of skin changes.

I am intended both for use by medical professionals (e.g., doctors, medical assistants) and for application in the cosmetic field.

I serve exclusively as a tool for skin analysis and am not responsible for making diagnoses.

Inappropriate use

If you use me in a manner other than described in this user manual, this does not correspond to my purpose. My manufacturer is not liable for damage resulting from improper use. In such cases, you, the user, bears the risk alone.

Even though I am protected against water and small particles, I am not waterproof.

Please do not use the attachment with glass (contact mode) on wounds; this is only permitted in contactless mode.

Description of the wireless connection

Wireless communication is carried out via the 2.4 GHz band, supporting the IEEE 802.11 b/g/n protocol.

My integrated module supports bandwidths of 20 MHz and 40 MHz. The transmission power is 20 dBm.

I have been designed to operate throughout the European Union without the need for special authorizations or usage restrictions.

You can use me freely in all Member States without having to worry about frequency approvals or special requirements.

Product overview

Accessories

- 1x USB-C cable for charging
- 1x attachment ring without glass for contact-free use

Materials

My elegant housing as well as the attachment rings are made of thermoplastic. The back plate and buttons are made of aluminum. My glass is high-quality quartz.

Inside me, you will find a circuit board, a lithium battery, a camera, and 16 concentrically arranged LEDs. Together, they ensure that I can deliver the best possible images of skin lesions.

Compatible devices

I am designed to operate independently of any specific end device. I transmit images wirelessly to your device.

You can connect me with any internet-enabled device, such as smartphones, PCs, or tablets. There are no restrictions regarding the operating system.

Minimum requirements for the terminal device

Wi-Fi connectivity: Your device must be able to connect to a Wi-Fi network in order to receive and transmit data from me.

Support of a modern internet browser: Your device must have an up-to-date web browser that supports modern web technologies.

Supported browsers (latest version required):

- Google Chrome
- Mozilla Firefox
- Microsoft Edge
- Apple Safari

Please ensure that your device complies with electromagnetic compatibility requirements to avoid interference.

Light modes and image effects

With the buttons on my back and the controls in my user interface, you can adjust the arrangement and color of my LEDs and apply different real-time effects to the image. This makes various skin details more visible.

Full circle

Homogeneous illumination of the skin area under examination

Half circle (left/right)

For detailed analysis of the skin surface

Warm white

Emphasizes deeper dermal tones

Useful for viewing underlying vascular structures but less effective for detecting superficial skin irregularities

Cool white

For general inspection and a natural rendering

Suitable for highlighting surface details and providing a clinical impression

Red

For examining deeper skin layers and assessing lesion depth

Enhances visualization of blood vessels, vascular lesions, and overall circulation

Less scattering, so surface details are less emphasized

Green

Assists in distinguishing between pigmentation changes and blood vessels

Provides strong contrast for superficial vessels

Suitable for highlighting small vessels and subtle erythema

Blue

For assessing the upper skin layers

Effective in identifying surface irregularities, scaling, and acne lesions

Skinverse® effect

Inverts light and dark areas to enhance contrast

Grayscale

Enhanced contrast, making differences in brightness and structure more apparent

Operation

Control & function check

Before each use, please perform a visual inspection for damage or wear and test if I function properly. You can recognize this by the fact that I connect to your device and the image quality does not deteriorate.

Do not use me if you detect any defects!

Charging the battery & application duration

I am equipped with a rechargeable lithium-polymer battery.

To charge me, connect me to a USB power supply using the supplied USB-C cable **4**. Use only CE-marked USB power supplies from reputable manufacturers that meet the safety requirements of IEC 60601-1 (medical electrical equipment).

My power consumption is 1–1.5 amps.

Please do not use power supplies or USB-C cables with visible damage.

My battery should be charged outside the patient environment (at least 1.5 m distance).

Before the first use, please fully charge me with the supplied USB-C cable.

Once the charging process is complete, my indicator light turns red.

I can be used continuously for up to 2 hours before I need to be recharged. You can see my current charge level at the top right of the user interface.

I can also be used while charging – however, this will extend the charging time.

Switching on & off

To switch me on, use the switch **1** on my side and push it upward (see illustration). You can tell I am switched on when my LEDs **2** light up in a circle in different colors.

Push the switch downward to switch me off.

I can start in two states:

- Normal start: I am ready for use and my LEDs light up fully in white.
- Low battery: I do not start and my LEDs flash red once every minute. This means I must be charged or switched off.

I have an integrated standby mode. If I am not moved for 30 seconds, I switch off automatically.

To switch me back on, simply move me briefly.

Connecting to your device

After switching me on, please go to the Wi-Fi settings of your multimedia device. My network "**Scaut**" will appear there.

Click on this network and enter the password "**doccheck**" for the first connection.

To connect your device with me, open a browser of your choice and enter the following URL: **http://scaut.local**


Once the connection is established, you can begin the examination.

Changing the password

When you use me for the first time, you will be prompted to change the Wi-Fi name and password to prevent access by others.

I will then restart. Reconnect to me by entering the credentials you saved.

If you ever forget your password – no problem! Press the top and bottom buttons simultaneously for 15 seconds to perform a reset. You can then reconnect with me using the login details mentioned above.

Button configuration

Capture: Press this button once to take a picture.

LED color: This button controls my LED colors. Each press cycles through the following colors:
Warm white > Cold white > Red > Green > Blue


LED array: This button controls my LED array configuration. Each press cycles through the following modes:

Full circle > Left half circle > Right half circle

User interface

After you have connected me to your multimedia device and entered the IP address in a browser, the following webpage will appear:

The drop-down menu **A** allows you to adjust settings such as changing the password or language, or performing an update. The following menu items are available:

- **Settings:** Change the SSID and password
- **Update:** Perform a software update (for more details, see the separate section on performing an update)
- **Languages:** Set the language of the user interface
- **About me:** Information about your Scout

Using icons **B**, you can change the shape and color of the scale. You can choose between a linear and a circular scale, each available in black or white.

With the magnifying glass **C**, you can zoom in on the area under examination to ensure that no detail is overlooked.

In the large window **D** in the center of the webpage, you will see the current live view.

Using the slider **E**, you can switch between different real-time image effects, such as Skinverse® and grayscale.

You can also access previously performed screenings. To do so, click on "Load Screening" **F** and select the desired file.

With the button **G**, you can start a new screening.

If the webpage becomes unresponsive, reload it with "Refresh" **H**.


At the top right, my current battery charge level **I** is displayed.

The buttons **J** on the right have the same functions as my physical buttons: taking a picture, adjusting the number of LEDs turned on, or changing their color.

Image capture

To take an image, you have two options: either press the capture button on my back or click the corresponding **J** button on your screen.

A preview of the captured image will appear with the following options:

If you want to properly document the image or take a series of images, click "Add to screening". This process is explained in the next section.

If you only want to save this single image, click "Save image". The image will be saved directly on your connected device with a generic name and the capture date.


If you do not want to save the image, click "Discard".

Note on image quality in live mode

Please note that image quality in live mode is lower than that of saved images. The final saved images have higher resolution and are therefore more suitable for detailed viewing and documentation.

Screening

If you started a "New screening" before capturing an image or clicked on "Add to screening" after capturing, the following pop-up window will appear:

Click on the body part in the illustration that matches the patient's body region from which you took the image.

You can change the gender of the figure using the corresponding icons. You can also rotate the figure in 90° steps, for example if the skin lesion is on the patient's back.

If you do not want to take any more pictures, click "Save screening".

If you want to capture more skin lesions, click "Next image". Take the image and click "Add to Screening" again. You only need to assign the body part.

When the examination is complete, click "Save Screening" – an HTML file named with the capture date will be

saved on your connected device. This file contains the patient data as well as the previously taken images. You can view them by clicking on the marked positions.

NOTE: The screen captures in this manual are for illustration purposes only. The actual display may vary depending on the software version.

Application on the patient

1. Positioning: Hold me securely by my handle and place me directly on the skin. Make sure I am centered over the skin lesion or the area to be examined. You can use me either with the contact glass or without it – depending on the type of examination you wish to perform.

NOTE: The camera focus is optimally adjusted and cannot be changed by turning the attachment ring.

2. Lighting: When you switch me on, I automatically start with cold white light. It is recommended to begin the examination with white light. Afterwards, you may select additional LED colors if needed to make further details visible. The exact configuration and sequence are described in the section above.

3. Capture: To take a picture, you can either press my top button or use the web interface. Ensure that I remain stable on the skin during image capture.

4. Abschluss: After the examination, carefully remove me from the skin. Instructions on cleaning and disinfection can be found in the chapter "Cleaning, disinfection and storage" of this manual.

Changing the glass attachment

I come with an attachment ring with glass pre-installed for contact mode. If you want to work in non-contact mode, unscrew the attachment ring (counterclockwise) and install the supplied ring without glass. Do not tighten the thread by force.

NOTE: The thread is not used for focusing; the attachment ring must be fully screwed on.

Changing the battery

1. Please remove my seven screws **3** using the allen key provided in the tuning kit.
2. Carefully turn me over and remove my back, including the buttons, from the rest of my body. Put it aside.
3. Next, remove the circuit board, including the battery and camera, from the housing.
4. Carefully disconnect my battery connector.
5. Now you can insert a new battery using the same principle. The connector only fits one way, so you can't go wrong.
6. Place the circuit board back into the designated position inside my housing.
7. Reassemble me by placing the aluminum housing with buttons back onto the rest of the housing and pressing it into place.
8. Finally, secure my housing again with the screws. Done – I'm ready for use again!

Updating the software

Software updates should be performed by you on the recommendation of my manufacturer. Before my first use, please check if an update is available.

To perform a software update, I need a Wi-Fi connection with internet access. I will automatically search for available Wi-Fi networks nearby. The easiest way is to use your smartphone's hotspot.

Please select your preferred Wi-Fi network and enter the corresponding password. I will then connect to this Wi-Fi and check for available updates, which will automatically download if present.

Please do not switch me off during the update.

If the update fails or is interrupted, restart me and repeat the steps above.

Software updates are regularly available – please check under the "Update" menu option.

Cleaning, disinfection and storage

Manual cleaning

Please clean me regularly. It's very simple: wipe me thoroughly with a damp cloth. My manufacturer recommends cleaning wipes with 70% isopropyl alcohol or soapy water.

If you use other cleaning or disinfecting agents, please test their compatibility with the materials first. Otherwise, my pleasant outer appearance may get damaged. Frequent cleaning doesn't bother me. I like to be clean.

Do not immerse me in cleaning solutions.

Sterilization

You CANNOT sterilize me. If you do, I will break and probably come to a sad end.

Storage

Store and use me only in a dry and dust-free environment. Please note the following limit values regarding temperature and humidity:

	In operation	Storage and transport
Temperature	0 °C to +35 °C	-20 °C to +70 °C
Humidity	45–80 %	10–90 %

After use, store me preferably in my packaging or in a suitable case.

Safety information

Usage restriction

I am not intended for invasive use.

I should also not be used in environments at risk of fire or explosion.

Please do not make any unauthorized modifications to me!

Light source warning

Do not look directly into the light source.

Do not use me for testing the pupillary light reflex, as this is not part of my intended purpose.

Data security

The protection of patient data is your sole responsibility (e.g., securing your smartphone with a PIN). Please ensure your device is adequately protected against third-party access.

All images are stored exclusively on your device.

You are responsible for creating regular backups to prevent data loss, for example in the event of device loss.

Notice to the user and/or patient

Any serious incidents related to the product must be reported to DocCheck Shop GmbH at lab@doccheck.com and the competent authority in your country.

Electromagnetic safety

I comply with the requirements for EMC emissions and immunity according to IEC 60601-1-2:2014 + A1:2020 and provide adequate protection for radio communication services.

My emission characteristics meet CISPR 11 Class A. I am suitable for use in industrial environments and hospitals.

Exclusion

I am not intended for operation in special environments as described in IEC 60601-1-2. These include in particular:

- Areas near MRI systems, which generate strong static magnetic and high-frequency fields.
- HF surgery environments where strong electromagnetic disturbances occur.
- Environments with intense electromagnetic fields, such as those produced by radiotherapy systems.

My emission characteristics make me suitable for industrial and hospital use (CISPR 11 Class A).

If I am used in a residential environment (for which CISPR 11 Class B is required), I may not provide adequate protection for radio communication services.

Precautions to avoid electromagnetic interference

To prevent undesirable events for patients or users caused by electromagnetic disturbances, please observe the following precautions:

- Use only accessories and cables specifically provided for me by my manufacturer. Other accessories may increase my electromagnetic emissions or reduce my immunity, leading to malfunctions.
- Portable devices using high-frequency communication (e.g., WLAN peripherals, antenna cables, or external antennas) should be kept at least 30 cm (12 inches) away from me.
- Avoid using me in close proximity to or stacked with other equipment, as this may cause improper operation. If such use is necessary, both I and the other equipment should be monitored to ensure proper function.
- Electrostatic discharges may cause interruptions in my lighting or the transmission of images to your device. Please use me in an environment that minimizes

such discharges, for example by ensuring appropriate flooring or humidity control.

Measures in the event of electromagnetic interference

If you notice that I cause or experience interference, you should take one or more of the following measures to resolve the issue:

- Reorient me or use me in another location.
- Increase the distance between me and the affected equipment.
- Connect me to a different power source while charging.

Essential performance characteristics

I do not have essential performance characteristics whose loss or impairment would present an unacceptable risk to patients, users, or third parties.

Compliance with emission and immunity standards

I comply with the requirements of IEC 60601-1-2 regarding:

- Emissions: CISPR 11 Class A, Group 1
- Immunity: see the following table:

Immunity testing	Basic EMC standard	EN/IEC 60601 Test conditions	Compliance
Magnetic fields at utility frequency (50/60 Hz)	EN/IEC 61000-4-8	30 A/m 50 and 60 Hz	30 A/m 50 and 60 Hz
Immunity of the enclosure port against nearby magnetic fields	IEC 61000-4-39	Test frequency: 30 kHz Modulation: CW Immunity test level: 8 A/m	8 A/m
		Test frequency: 134.2 kHz Modulation: Pulse 2.1 kHz Immunity test level: 65 A/m	65 A/m
		Test frequency: 13.56 kHz Modulation: Pulse 50 kHz Immunity test level: 7.5 A/m	7.5 A/m
Electrostatic discharge (ESD)	IEC 61000-4-2	Contact Discharge: ±8 kV Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Contact Discharge: ±8 kV Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV
High-frequency electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz–2.7 GHz 80% AM at 1 kHz	3 V/m 80 MHz–2.7 GHz

My immunity to high-frequency electromagnetic fields has been verified according to IEC 61000-4-3, Table 9.

Maintenance and warranty

I am maintenance-free.

If my glass or battery becomes defective, you can replace these components.

The warranty for the entire product will be voided if replacement parts are used that were not procured from the DocCheck Shop.

If you use me as intended and follow this instruction manual, the warranty for manufacturing and material defects is two years starting from the date of purchase.

Obvious abuse and accidental damage are excluded. Cosmetic damage, such as scratches on the surface, are also excluded from the warranty.

Disposal

Please note that I have to be sent to a separate collection point for electrical and electronic equipment.

If you have to dispose of the battery, this must be done at appropriate collection points. You should take a look at the specific disposal rules in your country for this.

Manufacturer

DocCheck Shop GmbH
Vogelsanger Str. 66, 50823 Cologne, Germany
doccheckshop.eu

Troubleshooting

Problem	Possible causes	Solutions
Scaut is not charging	- Defective USB-C cable or power adapter - Incompatible power adapter	- Test with another USB-C cable - Use a different power adapter
No connection to the terminal device or delayed image transmission	- Wrong network selected - Incorrect password - Distance too great - Interference from other devices	- Select the correct network - Check and re-enter the password - Restart the device - Perform a reset - Reduce the distance between me and your terminal device - Move away from possible sources of interference
No image visible / only gray screen	- Streaming page has frozen - Browser cache issue - Connection dropped	- Check Wi-Fi connection - Refresh the browser page - Use a different browser - Restart the device
LEDs do not light up	- I am not switched on - My battery is empty	- Switch me on - Check my battery status and recharge if necessary

Explanation of symbols

Medical device

CE marking

Manufacturer data (address)

Production date

Item number

Batch

Unique Device Identifier

Contains no latex/natural rubber

Permissible temperature range for storage and transport

Permissible humidity range for storage and transport

Fragile – handle with care!

Store in a dry place!

Direct current


Separate disposal of electrical and electronic appliances.

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